

MoGrip | Mobility Cane Attachment for Wrist Stabilization During Cane Use

04 | Team: Anna Betron, Ben Harnwell, Meera Rajasekaran, & Kittana Tran
Advisor: Dr. Sriram Balasubramanian, Associate Professor, Drexel University

Need

Mobility cane/white cane users experience mild/moderate **wrist pain** during cane use. There are no existing solutions that **prevent wrist discomfort** without **compromising functionality**.



*Ambutech No-Jab Mobility Cane


Objective

Our objective is to design a **secure, portable cane attachment** to **reduce wrist adduction** and **reduce wrist flexion/extension** to minimize pain during cane use.

Design Inputs



- Encourage a **neutral wrist position**, reducing **adduction & flexion/extension** while maintaining sweeping motion
- Fits within **portable** profile of cane
- Attaches to existing cane **securely**

Solution




Intended Use

Muscle Groups Utilized During Use




Verification Testing

Wrist Angle Testing - PASSED
Wrist adduction, extension & flexion were reduced while maintaining sweeping distance



Secure Attachment - PASSED

Achieves a fully compressed No-Jab spring with & without the MoGrip



Portability Testing - PASSED

Device is within the correct dimensions to be portable

Dimension	Criteria	Results
Mass	< 81.83 g	72 g ✓
Handle Length	6.5 to 10.5 cm	9.7 cm ✓
Head Width	< 6.25 cm	6 cm ✓

Conclusion

- **Fully 3D Printed/Printable:** \$0.92 /device
- Comfortable wrist position that **reduces pain**
- Improves **quality of life** and **user experience**

Future Work

- Explore **other materials and manufacturing methods**
- Improve **strap clasp** design
- Incorporate a not-in-use **carrying method**

Scan QR Code to Watch a MoGrip Demo Video!

